

# OrganoStation

Welcome electrophysiology in brain organoids



Precise electrode location



Fluorescence imaging



Optogenetic stimulation



Acute drug application



# Technical specifications

- ✓ Simultaneous acquisition of up to 512 electrical signals at 30 kHz
- ✓ Micromanipulators to allocate the electrodes with micrometer precision
- ✓ Long working distance infra-red and fluorescence wide-field imaging
- ✓ Electromagnetic isolation and chemically resistant components
- ✓ Peristaltic drug application and temperature control of sample







# Precise electrode location

Precise movement of neural probes together with tip sharpening permits penetrating the organoid with minimal damage.


The microelectrodes are arranged in **200-300  $\mu\text{m}$**  shanks to sample the outer core of the brain organoid – where active neurons exists.




## Application example

*Low-impedance small microelectrodes are designed to record single-active neurons or neuronal population events in brain organoids.*

### Single-neuron activity




### Population event





# Precise electrode location

- ✓ Precise ( $< 1 \mu\text{m}$ ) motorized (4.8 nm resolution) XYZ micromanipulators with 25 mm travel distance.
- ✓ XYZ manual controlled (25 mm) stage oxidized to reduce toxicity or chemical degradation.
- ✓ Vibration-damping optical table (600 x 900 mm L x W) and the associated passive-isolation table.








# Optogenetic stimulation

The heterologous expression of an optogenetic actuator (a channelrhodopsin, e.g. ChR2) permit using light to activate neurons in brain organoids and test its function.





# Fluorescence imaging

- ✓ Large (20-80 mm) working distance and field of view (12 x 9 mm) with 0.6-1.7X magnification range.
- ✓ Computer vision infrared video-camera with 3.5  $\mu\text{m}$  per pixel resolution up to 227 frames per second.
- ✓ Broadband or single-colored (e.g., 470 nm) strong up-right illumination (up to 1200 mW).
- ✓ Filter sets for green (470 BP/495 LP/525 BP) or red (590 BP/660 LP/700 BP) fluorophores.







# Acute Drug Application


Precise **temperature control** for solution and organoid to preserve physiological conditions.

Recirculation of small volume for **cost-effective** drug application.



## Application example

*Effect of a glutamate antagonist (CNQX) and GABA antagonist (Gabazine) on extracellular spikes recorded in human brain organoids.*





VIENNA SCIENCE & ENGINEERING

[www.vise-og.com](http://www.vise-og.com)